**Software UniChrom to unify work
with different analytical equipments
from different vendors**

There is the following situation in the researches laboratories: very heterogeneous park of analytical equipments

From different Vendors

Gas chromatographs

Liquid chromatographs

**Spectrometers
UV, IR, Vis, AES, AAS, MS**

Device model	Vendor	Software name
Agilent (HP) -USA	Agilent (Hewlett-Packard)	ChemStation
Varian - USA	Varian	Star WorkStation / Galaxy
Perkin-Elmer - USA	Perkin-Elmer	TurboChrom / TotalChrom
Shimadzu GC17/2010 -Japan	Shimadzu	GC Solution / Class VP
Crystall 2000M/5000 -Russia	JS «Chromatek»	Chromatec-Analytic
Crystalluks 4000 Russia	JS «Metachrom»	NetChrom
Tswett-800	Tswett Ltd.	Tswett-Chrom
MiLiChrom-5	JS «Nauchpribor»	E-Chrom
MiLiChrom A-02	JS «EcoNova»	M-Chrom
ThermoFinnigan Trace2000/Focus	Scientific Software	EZChrom / Xcalibur

The analogue situation was about 30 years ago.

There were many text processor packages in Laboratories.

Demand to free migration of text documents between different computers has generated the following situation:

- the computers may be different but the text processor software Microsoft Word is the same in the most offices computers over the world.

All devices are controlled in unified style

UniChrom™ - driver installer and configuration editor


File Edit

For driver installation just drag with the mouse Instrument icon from the "Available section" (at the right) to "Installed" section (at the left). For driver Removal drag the Instrument icon from the "Installed devices" section to the "Available" instruments. Instrument settings is done in Settings Dialog on double click on Instrument icon in left section.

Installed devices	Available for installation
 Agilent 7890.Agile...  Cluster.Cluster	 Acme-6000 (acm6k.inf)  Agilent 7890 (ag7890.inf)  Aquilon ADC (AquilADC.inf)  Chromat-900 (chromat9.inf)  Cluster (hplcc.inf)  Counter (counter.inf)  E-24 (e24.inf)  EMCES (emces.inf)  Fluorat-02-02M (flu02.inf)
 GC20.GC20  Trace-2000...	 GC17A-2 (gc17a2.inf)  GC20 (gc20.inf)  GC-Sim (gcstub.inf)  HP5890 (hp5890.inf)  HP6890N (hp6890.inf)  K-2xxx (k2.inf)  Kapel (kapel.inf)  Lambda-1010 (lambda.inf)  LNet (lnet.inf)
 Varian cp3800.Vari...  Кристалл-5...	 LNetI (LNetI.inf)  LNetI-800+ (lnet800.inf)  LNet-X new (lnetx.inf)  LZI (lzi.inf)  MLC-16 (mlc16.inf)  MSS-ADC (mssadc.inf)  OxyText (oxy.inf)  PeriChrom PR2100 ...  PLC-24 (plc.inf)
	 Polaron-32 (polaron.inf)  Riftek (riftek.inf)  Solar AP-2110 (ap2110.inf)  SoundBlaster (sbwin32.inf)  Trace-2000 (trace2k.inf)  Varian cp3800 (cp3800.inf)  Z-Lab (zlab.inf)  АПК-02 (adc2.inf)  Кристалл-2... (cr2000m.inf)
	 Кристалл-5... (cr5000.inf)  Кристалл-P... (cr5000p.inf)  КристалЛю... (cr4000.inf)  ЛГХ-3000 (lgs3000.inf)  ЛХМ-2000 (lcm2000.inf)  Марафон-II (marathon2.inf)  МилиХром (milichrom.inf)  МиЛиХром... (mlca02.inf)  Стайер (marathon.inf)
	 Стайер-Grad2 (stayerg.inf)  Стайер-ISO (stayeri.inf)  СтайерКД (stayercd.inf)  Хромос  ГХ-1000 [...]  Цвет-800+ (c800.inf)  Яуса (yausa.inf)

Different UV, IR, AAS and AES

UniChrom – AES

Distinctives of Unichrom

The screenshot displays the NAS UniChrom V software interface, which is used for gas chromatography data analysis. The interface is divided into several sections:

- Instrument Parameters:** A central panel shows various instrument settings for the GC Instrument. Parameters include Oven Temp (50.0), Inj1 Carr (0.0), Inj1 Split (20.0), Inj2 Carr (0.0), Inj2 Split (20.0), Det2 Temp (240.0), PID Temp (150.0), PID Gain (0.0), and Events List (0.0).
- Temperature Profile:** A graph on the right shows the oven temperature profile over time, starting at 50.0, ramping up to 200.0, and then holding constant.
- Method Configuration:** A section titled "Object: Chromatograph / Oven / Temperature" shows the chosen mode: "4 - ГОСТ 23781 - 5 изо - прорг". It includes a tree view of the method structure and a table of parameters.
- Method Table:** A table titled "Programme: 40 min" lists the temperature program steps:

Rate	Value	Interval
1	0	5
2	10	220
3		18
- Chromatogram:** A large plot at the bottom shows the chromatogram with peaks labeled with chemical names and their retention times. The x-axis is labeled "s" and ranges from 0 to 2000. The y-axis is labeled "0.2".
- Library:** A small window on the right shows a library search result for "Fragm.dY" with a value of "0,7093".
- Navigation and Tools:** The top of the interface features a menu bar (File, Edit, View, Tools, Windows, Help) and a toolbar with various icons for file operations and analysis.

Abstract representation of the method focused the Analyst to analytical details instead of instrument peculiarities

Quality control of **Spirit Drinks**. Thousands food control laboratories: over the world, day-and-night, there are many different analytical equipments but the software may be the same.

UniChrom - проверка системы - Windows Internet Explorer

http://unichrom.com/vodka/ethanol-solutionr.shtml

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Free migration of methods between different instruments.

Unichrom Project Road Map

1. Certification of the software package UniChrom.
2. Development and Approval of the unified State Standard for validation different analytical equipment with UniChrom.
3. Introducing in the Accredited testing Laboratories
4. The biggest Oil Refineries
5. Food control Laboratories
6. Toxic laboratories in narco dispensaries
7. Studies and Researches Laboratories in the leading Universities
8. Participation in leading international conferences and exhibitions.
9. Ready analytical solutions and there openness
10. Openness presentation in *Internet*.

Control laboratory in JS «Nevynnomusky Azot» (Russian Federation) . **UniChrom** from one PC controls simultaneous: **Crystal-5000** (3), **Crystall-2000M** (3), **HP6890N** (1), **Tswett-800** (2), **Tswett-100** (3).

Day-and-night Unichrom is working in the following biggest petrochemical enterprises:

- **JS “Novopolotsk Oil Refinery”**,
- **JS “Mosyrskii Oil Refinery”**,
- **JS “Lisichanskii Oil Refinery”**,
- **JS “ArcelorMittal Kryviy Rih”**,
- **JS “Mogilov HIMVOLOKNO”**,
- **JS “Grodno Azot”**,
- **JS “Nevynnomuskii Azot”**,
- **JS “Schekino Azot”**,
- **JS “Cherkaskii Azot”**,
- **JS “Severodonetskii Azot”**,
- **JS “Polimir»**,
- **FSI “Siberian Chemical Enterprise” (Russia) and Busher NC (Iran).**

List of the biggest petrochemical enterprises

UniChrom (tm) GC instrument support matrix - Mozilla Firefox

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Categories / Chromatography / UniChrom (tm) GC instrument support matrix

Solutions: [Gas Chromatography](#) [Liquid Chromatography](#) [Capillary Electrophoresis](#) [ADC for Chromatography](#) [Addition](#)

Applications fields: [UniChrom](#) [Restek Products](#) [Product Support](#) [Download](#)

The unified approach to GC instrument control makes the understanding your system easy and clear. Graphical representation of system components like detectors and detectors either as chromatography aspects, the optimal chromatography conditions make the system easy to control and fit your needs.

Set of different GC instrument parameters is collected in a single "Mode" and each method can easily employ different mode for every chromatography in batch.

UNICHROM - SUPPORTED GAS CHROMATOGRAPHIC INSTRUMENTS

Agilent 7890a.
 Producer: Agilent Technologies. UniChrom supports 7890 instruments with automated liquid samplers. Simply the best in the world.

Agilent/HP-6890, 6890 Plus, 6890N, Agilent 6890.
 Producer: Agilent Technologies. UniChrom supports all the family of 6890, 6890 instruments with 7673 and 7683 automated liquid samplers. Simply the best in the world.

Agilent/HP-6890, HP-6890A, HP-6890 Series II, HP-4890D.
 Producer: Agilent Technologies, Hewlett-Packard. Well known GC instruments which often is mentioned as key tone of computer chromatography. Still works in many laboratories.

Shimadzu GC-2010.
 Producer: Shimadzu Inc. Honorable response to the competitor, Baazal. Up to 9 atmosphere at column inlet. Japanese quality.

Shimadzu GC-17A.
 Producer: Shimadzu Inc. Simple and robust workhorse for laboratory. Since instrument does not have its own built-in ADC, the choice of it is less attractive. We support these instruments thru our ADC LIBNET, which gives the GC-17A wide range measuring capability for ionization detector.

Crystal-2000.1
 Producer: JSC Chromatic. The new generation chromatograph. Simply the best in Russian market. With its characteristics comparison with PE Clarus, AutoSystem XL, Thermo FluoTag. [Trace 2000](#). Right price / performance ratio.

Crystal-2000.2
 Producer: JSC Chromatic. Compact gas chromatograph built upon Crystal-2000.1 platform with the set of improvements and new features.

Crystal-2000M.
 Producer: JSC Chromatic. Gas chromatograph with automated liquid sampler. Suitable for most types of analysis. Full UniChrom control. Simply the best in Russia. Now it is preceded by newer model [Crystal-2000M](#).

Varian cp-3800
 Producer: Varian Inc. Full instrument control except configuration editing (easy done from instrument keyboard).

Trace 2000
 Producer: Thermo FluoTag Inc. Among the top GC of the world but [Crystal-2000](#) is better. Electronic flow pressure controller (DPFC - dynamic P/F controller).


Crystal-4000
 Producer: JSC Chromatic and JSC (Meta-Chrom). Simple chromatograph with electronic temperature and flow controller. UniChrom driver supports full control of this instrument, and also ramped programs for all temperature zones and all flow pressure controllers. The only advantage - its price. Its tasks can be compensated by personal qualification.

Any gas chromatograph with analogue output or **directly** from detector.
 Virtually any GC instrument can be connected to computer with the help of ADC [LIBNET](#).

GC-Sim. Gas chromatograph simulator. Just training driver. Download UniChrom, install GC-Sim and enjoy the learning process.


Download Center is opened Day-and-Night

Unified style for different operation systems

Windows Vista 32bit

Debian GNU/Linux

Results:

Certification in Belarus, Russia and Ukraine

Day-and-Night service at the biggest petrochemical enterprises:

JS “Novopolotsk Oil Refinery”, J “Mosyrskii Oil Refinery”, JS “Lisichanskii Oil Refinery”, JS “ArcelorMittal Kryviy Rih”, JS “Mogilov HIMVOLOKNO”, JS “Grodno Azot”, JS “Nevynnomuskii Azot”, JS “Shchekino Azot”, JS “Severodonetskii Azot”, JS “Polimir”, Federal State Int “Siberian Chemical Enterprise” (Russia) and Busher NC (Iran).

More than 1 million **alcohol in blood analysis (very social important)** have been done using Unichrom in Belarus and in the city of Moscow.

There is Device Development Kit (DDK) on the site www.unichrom.com. So independent researches can develop drivers for new analytical devices independently.

Software package *UniChrom* is proposed for High Schools on free charge basis.

Plans:

To expand park of new analytical devices under control of Unichrom.

Total support of modern multidimensional detectors (UV/Vis DAD, GC/MS, GC x GC)

To supply modern Laboratory Information Management Systems (LIMS) with genetic plane texts **E-Lab**

Installation with local market language: Chinese or Korean

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